

INTRODUCTION

1. AutoZyme Creatinine is a reagent set for determination of creatinine based on **initial rate method** using Alkaline Picrate.
2. AutoZyme Creatinine is a **single reagent system**, using one step procedure.
3. AutoZyme Creatinine has **one step reconstitution**. It involves mixing of Picrate and Diluent reagent.
4. AutoZyme Creatinine is a **High Stability Reagent**.
5. AutoZyme Creatinine is **linear** upto 30 mg%.
6. Creatinine can be determined in **180 seconds**.
7. AutoZyme Creatinine can be used on any **Spectrophotometer, Discrete semiautomated and Automated analyzer**. Programme can be designed for any specific analyzer upon request.

PRINCIPLE

Creatinine in alkaline medium reacts with picrate to produce orange colour. This colour absorbs light at 492 nm.(490 - 510 nm). The rate of increase in absorbance is directly proportional to the concentration of creatinine in specimen.



PREPARATION OF WORKING SOLUTION

Prepare working solution by mixing equal volume of Picrate Reagent and Diluent Reagent.

REAGENT STORAGE & STABILITY

The reagents are stable till the expiry date stated on the bottle label, when stored at R.T. (25-30°C).

The working solution is stable for 30 days at 2-8°C.

COMPONENTS & CONCENTRATION OF WORKING SOLUTION

Component	Concentration
• Sodium Picrate	7.7 mmol/l
• Sodium Hydroxide	500 mmol/l

SPECIMEN COLLECTION & PRESERVATION

Blood should be collected in a clean dry container. Avoid use of plastic or siliconized container which may prolong clotting time. Samples should not be collected during PSP/BSP clearance test. For plasma separation Heparin (200 IU/ml blood) may be used as anticoagulant.

Creatinine in serum and plasma is stable for 2 days when stored at 2-8°C.

PROCEDURE

- Reaction type Initial rate
- Reaction direction Up
- Wavelength 492 nm.(490 - 510 nm.)
- Flowcell temperature 30°C / 37°C
- Zero setting with Distilled water
- Delay time 30 seconds
- No. of readings 2
- Interval 60 seconds
- Sample volume 0.05 ml (50 µl)
- Reagent volume 1.0 ml
- Standard concentration 2 mg %
- Factor $2 \div \Delta \text{ Abs. of standard}$
- Linearity 30 mg/dl

Manual assay procedure

Prewarm the required amount of working solution to 30°C/37°C before use.

1.0 ml procedure

Standard / Sample	0.05 ml (50 µl)
Working Solution	1.0 ml

Mix and start stopwatch simultaneously. Record absorbance assay mixture at exactly 30 seconds after Standard / Specimen addition and then again at 90 seconds.

Note : It is recommended to run a creatinine standard with each batch of assay.

Calculation:

Calculate the average change in absorbance per minute ($\Delta \text{ Abs.}$) of standard & specimen(s).

$$\Delta \text{ Abs.} = \text{Abs. at 90 sec.} - \text{Abs. at 30 sec.}$$

$$\text{Serum Creatinine (mg\%)} = \frac{\Delta \text{ Abs. of Specimen}}{\Delta \text{ Abs. of Standard}} \times 2$$

