


INTRODUCTION

1. AutoZyme **STAT** Glucose is a reagent set for determination of **True Glucose based on enzymatic method** using Glucose oxidase and Peroxidase.
2. AutoZyme **STAT** Glucose estimates glucose in just **one minute** using an initial rate method and in 7 minutes by end-point method.
3. AutoZyme **STAT** Glucose pack sizes 10 x 100 ml and 4 x 500 ml are provided with **2 standards** of 100 mg% and 300 mg%.
4. AutoZyme **STAT** Glucose is **Linear** upto 700 mg% glucose for kinetic assay procedure and 500 mg% using end-point procedure.
5. AutoZyme **STAT** Glucose is a High Stability Reagent.
6. AutoZyme **STAT** Glucose can be used on any **Spectrophotometer, Discrete semiautomated and Automated analyzer**. Programme can be designed for any specific analyzer upon request.
7. AutoZyme **STAT** Glucose has **one step reconstitution**. It involves the mixing of ENZYME and DILUENT.
8. **Sodium Fluoride** (as an anticoagulant upto 10 mg/ml blood) does not effect glucose assay.
9. The influence of **Ascorbate, Billrubln, Antidiabetic drugs and Haemoglobin** is negligible.
10. **Very small volume** of serum or plasma is required for assay.
11. The **shelf-life** of AutoZyme **STAT** Glucose is 18 months.

PRINCIPLE

Glucose oxidase (GOD) converts Glucose to gluconic acid. Hydrogen peroxide formed in this reaction, in presence of Peroxidase (POD), oxidatively couples with 4-aminoantipyrine / phenol to produce red quinoneimine dye. This dye has absorbance maximum at 505 nm. (500-550 nm.). The intensity of the colour complex is directly proportional to the Glucose in specimen.



REAGENT AVAILABILITY

	10 x 100 ml (Code : GU-2)	4 x 500 ml (Code : GU-3)
• ENZYME Vials	10	4
• STANDARD bottle - 100 mg%	1	1
• STANDARD bottle - 300 mg%	1	1
• WORKING Solution Container*	1	—
• DILUENT bottles	2	4

*(Empty bottle for preparation and storage of working solution)

*GLUCOSE standard 100 mg% and 300 mg% are assayed against National Bureau of Standards (NBS) reference material from Washington.

ENZYME and STANDARD are to be stored at 2-8°C.

Diluent Reagent may be stored below 30°C and away from direct light.

The reagents are for *in vitro* diagnostic use.

PREPARATION OF WORKING SOLUTION

4 x 500 ml (Code GU-3)

Transfer the content of one vial of **Enzyme** into one bottle of **Diluent** and mix by gentle swirling or inversion. Do not shake vigorously.

10 x 100 ml (Code GU-2)

Transfer 100 mL. **Diluent** into working solution bottle. Now add contents of one vial of **Enzyme** into it. Mix by gentle swirling or inversion. Do not shake vigorously.

REAGENT STORAGE & STABILITY

Stability of reconstituted working reagent is printed on the bottle label. **(DO NOT FREEZE)**. The working solution should be stored in the dark bottle (diluent bottle or the working solution bottle) provided. This is critical because the reagent is light sensitive (auto oxidation of chromogen system by light and air), it should therefore be kept away from direct light.

COMPONENTS & CONCENTRATION OF WORKING SOLUTION

Component	Concentration
• Phosphate Buffer; pH 7.0	170 mmol/l
• Glucose oxidase	15000 IU/l
• Peroxidase	1500 IU/l
• 4-aminoantipyrine	0.28 mmol/l
• Phenol	16 mmol/l
• Stabilizers and inactive ingredients	

SPECIMEN COLLECTION & PRESERVATION

Blood should be collected in a clean dry container. Serum or plasma should be separated from the cells at the earliest possible (within 30 minutes), as the rate of glycolysis is approximately 7 mg% per hour at room temperature.

For plasma separation following anticoagulants may be used.

• EDTA	2 mg/ml of blood
• CITRATE	6 mg/ml of blood
• HEPARIN	200 IU/ml of blood
• OXALATE	3 mg/ml of blood
• SODIUM FLUORIDE	10 mg/ml of blood

Sodium Fluoride is preferred as anticoagulant due to its antiglycolytic activity. The higher concentration of Sodium fluoride i.e. more than 10 mg/ml blood should be avoided as it may inhibit the colour development. Glucose is stable for 24 hours in neatly separated plasma and serum. If the estimation is not possible within 24 hours then the specimen should be preserved at -10°C and should be used within 30 days.

PROCEDURE FOR END-POINT

- Reaction type.....End-Point
- Reaction time7 mins. at 37°C/15 mins. at R.T. (25-30°C)
- Wavelength505 nm. (500-550 nm.)
- Zero setting with.....Reagent Blank
- Blank absorbance limit < 0.300 Abs.
- Sample volume.....0.01 ml (10µl)
- Reagent volume.....1.0 ml
- Standard concentration100 mg%
- Linearity500 mg/dl

Manual assay procedure

Prewarm at room temperature the required amount of working solution before use.

Perform the assay as given below :

1.0 ml procedure

	Serum / Plasma	Standard	Blank
	0.01 ml	0.01 ml	—
Working solution	1.0 ml	1.0 ml	1.0 ml

Incubation

Incubate the assay mixture for 7 minutes at 37°C or 15 minutes at room temperature (25 - 30°C). After completion of Incubation period measure the absorbance against blank at 505 nm. Final colour is stable for at least two hours if not exposed to direct light.

Calculation :

$$\text{GLUCOSE In mg\%} = \frac{\text{Absorbance of Sample}}{\text{Absorbance of Standard}} \times 100$$

